

HOSPITAL CLOTH



ITEM	TEST DATA	
Fibre Composition:	50% viscose	
	50% polyester	
Bonding Method:	Needle punch on fibre web	
Mass per Unit Area:	45g/m ²	
Breaking strength, N, min machine direction (MD), dry, wet cross direction (CD), dry, wet:		
	Dry MD breaking strength \geq 90 (N/5CM)	111.7 N/5CM
	Dry MD breaking strength \geq 75 (N/5CM)	98.5 N/5CM
	Wet MD breaking strength \geq 75 (N/5CM)	80.7 N/5CM
	Wet MD breaking strength \geq 90 (N/5CM)	65.2 N/5CM
Absorbency:	Absorbent \geq 400 (%)	Dry weight: 4.77g
		Wet weight: 22.35g
		469% > 400%
Colours:	Blue, Green, Red, Yellow, White	